


DATA SHEET

Binder Recovery Apparatus by Vacuum



75-B0025/B

General description

Used to remove the solvent from the binder/solvent solution in order to determine directly the total content binder in the aggregate/binder mixtures.

The apparatus consists of a power operated vacuum pump, fit with vacuum regulator, producing a vacuum down to 200 mbar, a thermostatically controlled water bath, and two flat-bottomed flasks 250 ml capacity with rubber bungs and connections, All necessary fittings and connections complete the set.

The water bath can also be used for other applications as for example for the Marshall test. For more information see Water baths, model 76-B0066/A. For more information on the Vacuum pump, see Vacuum pumps, model 86-D2001 with vacuum gauge regulator 86-D2004/5.

- Power rating: 1380 W (1200 water bath, 180 vacuum pump)
- Weight approx.: 23 kg

Standards

- EN 12697-1

Products

75-B0025/B

Binder recovery apparatus, vacuum pump method. 230V/50-60Hz/1ph

75-B0025/BZ

Binder recovery apparatus, vacuum pump method. 110V/60Hz/1ph